

SeaMed[®]



LINEAR

WIRELESS POCKET ULTRASOUND

CE FDA ISO

 DENİZLER
MEDİKAL

More Professional Imaging with
192 Element Probe



Heat-resistant Special
Cooling Coating

Non-Slip Grip Pads


Longer Battery Life with
4200 mAh Dual Battery Cells



Compatible with Android, iOS, and Windows operating systems.
Can be used on PC, tablets, and smartphones.

 IP68

IP68 level water and dust resistance. The body provides maximum protection against dust and impacts, designed to protect the device even during prolonged exposure underwater.



TECHNICAL SPECIFICATIONS

Probe: Wireless Linear Probe

Elements: 192

Frequency: 7.5-10Mhz

Depth: 20-80mm

Presets: Thyroid, Small parts, Pediatrics, Vascular, Carotid, Breast, MSK, Nerve

Host: IOS/Android/Windows, Tablet, Smart phone, PC

Connection: WiFi

Software Adjustable Parameters: GN(gain), D(Depth), ENH(Enhancement), DR(Dynamic range), F(Frequency), FocusPos, PRF, WF, Mode, 8TGC, Biopsy, Annote

Measurement

B: Length, Area, Circumference, Angle, Trace, Distance

B+M: Heart rate, Time, Distance

B+PW: Velocity, Heart rate(2), S/D, Depth

Playback Frames: 100, 200, 500, 1000 optional

Display Mode: B, B/M, Color, PW, PDI, B+B

WiFi Type: Built-in WiFi, 802.11g/20MHz/5G/450Mbps

Channel: 64

Frame: 24s/f

Gain: 40-110 Db

Power: by built-in 4200mAh battery

Battery Replaceable: Yes

Battery Last: 3hours (Working time), 12hours (Stand by time)

Charger: Wireless charger, USB cable

Cooling Coating: Yes

Waterproof: IP68

Anti-skid design: Yes, 4 sides